

## Affordable hand protection, offering excellent comfort

- **Comfortable fit:** This disposable, nitrile glove has a thin formulation for reliable protection, plus improved comfort and flexibility
- **Secure grip:** The textured fingertips of the glove provide an improved grip for a variety of applications
- **Smart packaging:** Packed with 300 gloves per dispenser box\* to provide additional value and convenience

\*250 pieces/dispenser for size XL



### Industries

- Food Processing
- Automotive
- Life Sciences

### Recommended For

- Catering/Meal preparation
- Food processing and handling
- Inspection of parts, equipment
- Assembly and inspection of components
- Equipment repair and maintenance
- Picking and fastening components
- Adjusting parts and systems
- Handling of paint/glue guns and sprayers
- Transferring liquids and solids
- Sample taking and processing

TECHNICAL DATA SHEET

|                        | Product Information   |
|------------------------|---|
| External Glove Surface | Textured Fingers  |
| Audit Standards        | ISO 13485, ISO 9001   |
| Product Certification  | Medical Device Regulation (EU) 2017/745 Class I, Personal Protective Equipment Regulation (EU) 2016/425 Category III risks          |
| Packaging Overview     | 300 gloves per dispenser, 10 dispensers per shipper, 3000 gloves per shipper<br>**Size XL 250 pieces/dispenser, 2500 pieces/shipper |
| Storage Instructions   | Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.                            |
| Country Of Origin      | China   |
| Product Segmentation   | Standard Thin   |
| Antistatic             | Not Tested  |


PHYSICAL PROPERTIES

|   | Typical Values |             | Testing Method                     |
|---|----------------|-------------|------------------------------------|
| Length (mm/inches)                      | 230 / 9.1      |             | ASTM D3767,EN 420                  |
| Freedom from Holes (Inspection level I) | 1.5 AQL        |             | ASTM D6319-10,ASTM D5151-06 (2011) |
| Palm Thickness (mm/mils)                | 0.06/2.4       |             | ASTM D3767,EN 420                  |
| Finger Thickness (mm/mils)              | 0.08/3.1       |             | ASTM D3767,EN 420                  |
|   | BEFORE AGING   | AFTER AGING |                                    |
| Ultimate Tensile Strength (MPa)         | ≥ 14           | ≥ 14        | ASTM D6319-10                      |
| Elongation at Break (%)                 | ≥ 500          | ≥ 400       | ASTM D6319-10                      |
| Force at break (N)                      | 3              | 3           | EN 455-2                           |

ORDERING INFORMATION

| Size         | XS (5.5 - 6) | S (6.5 - 7) | M (7.5 - 8) | L (8.5 - 9) | XL (9.5 - 10) |
|--------------|--------------|-------------|-------------|-------------|---------------|
| Product Code | 82133060     | 82133070    | 82133080    | 82133090    | 82133100      |

Performance Standards and Regulatory Compliance



For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

|  |   |  |
|--|---|--|
| <b>Europe, Middle East &amp; Africa Region</b><br>Ansell Healthcare Europe NV<br>T: +32 (0) 2 528 74 00<br>F: +32 (0) 2 528 74 01<br><br><b>Asia Pacific Region</b><br>MY +603 8310 6688<br>apac.medical@ansell.com<br>CN: +86 21 38275000<br>infochina@ansell.com<br>JP: +813 5549 8151<br>info.medical.jp@ansell.com | <b>North America Region</b><br>Ansell Healthcare Products LLC<br>T: +1 800 800 0444<br>F: +1 800 800 0445<br><br><b>Latin America &amp; Caribbean Region</b><br>BR +55 (11) 3356-3100<br>latam.medical@ansell.com<br>MX: +52 442 296 2050<br>latam.medical@ansell.com<br>CO: +57 1 288 3247<br>latam.medical@ansell.com | <b>Australia</b><br>Ansell Limited<br>T: +61 1800 337 041<br>F: +61 1800 803 578<br><br><b>Russia</b><br>Ansell PYC<br>Ten. +7 495 258 13 16 |
|--|---|--|

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending: [www.ansell.com/patentmarking](http://www.ansell.com/patentmarking) © 2024 Ansell Limited. All Rights Reserved.

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

